

# LATEX SURGICAL GLOVES (STERILE) GUANTES PARA CIRUJANO DE LATEX ESTÉRILES



## Latex Surgical Gloves (Sterile)

KNOXVILLE MEDICAL Latex Surgical Gloves are traditionally known for excellent protection, fit and touch sensitivity. They are manufactured with naturally curved, tapered fingers which provide comfort and reduce hand fatigue. KNOXVILLE MEDICAL Surgical Gloves also feature an inverted bead cuff that reduces the chance of contamination due to roll-down. The unique package for KNOXVILLE MEDICAL Surgical Gloves is guaranteed not to tear, yet opens easily and cleanly. This package is impervious to contamination and minimizes latex deterioration due to ozone aging. KNOXVILLE MEDICAL Latex Surgical Gloves are visually inspected and meet or exceed government quality standards.

- Excellent comfort and touch sensitivity
- Curved, tapered fingers for longer wearing or hard-to-fit hands
- Longer cuff and ring roll design add extra protection
- Conveniently packaged in sterile pairs

## Guantes para Cirujano de Latex Estériles


Los guantes de cirujano estériles KNOXVILLE, son tradicionalmente conocidos por su excelente protección, suavidad y su adaptación a la mano. Son fabricados con la curva natural de la mano, dedos de fácil manipulación que provee confort y reduce la fatiga de la mano.

Los guantes de cirujano estériles KNOXVILLE, introducen un aro de latex que reduce la posibilidad de contaminación durante el proceso de entrada a la mano.

El empaque de nuestros guantes está garantizado a no romperse fácilmente, fácil de abrir. Este paquete está garantizado contra contaminación y minimiza el deterioro del latex por el O<sub>3</sub>.

Los guantes de cirujanos estériles KNOXVILLE son visualmente inspeccionados y cumplen con las normas de calidad standards internacionales.

- Excelente confort y suave sensibilidad.
- Dedos curvados, de fácil manipulación para largos períodos de uso.
- Diseño de manga larga y aro de latex, agregan una extra protección.
- Empacados en juegos de 1 (par) estériles.



### ORDER FORM / FORMA PARA ORDENARLOS

KVG 65	size 6-1/2	24 pairs/box; 15 boxes/case	360 und
KVG 70	size 7	24 pairs/box; 15 boxes/case	360 und
KVG 75	size 7-1/2	24 pairs/box; 15 boxes/case	360 und
KVG 80	size 8	24 pairs/box; 15 boxes/case	360 und
KVG 85	size 8-1/2	24 pairs/box; 15 boxes/case	360 und